

Sensigent eNose[®] Sensors

Medical Applications of Nanocomposite Sensor Arrays



Product Platform

Cyranose[®] 320




- Handheld Chemical Vapor Detector
- Designed for Industrial / Commercial Use
- Evaluated as an exempt Research Tool to determine efficacy of sensor arrays for a variety of medical applications

Nanocomposite Sensor Technology



Vapor passes over the polymer matrix and produces a change in dc resistance for each sensor

32 chemical sensors in standard array




Using pattern matching algorithms, the data is converted into a unique response pattern (PCA, CDA, ANN, SVM)

Medical Applications

**Premise –
Routine Diagnosis
and
Health Monitoring
through Breath Analysis**

Low cost and low power nanotechnology sensors will enable affordable and reliable devices for home health and point of care products

**product analogy:
spirometer on a PC card**



History of Breath Analysis

Hippocrates - treatise on breath aroma and disease

Lavoisier and Laplace (1784) - showed that respiration consumes oxygen and eliminates carbon dioxide

Nebelthau (mid 1800s) - showed that diabetics emit breath acetone

Anstie (1874) - isolated ethanol from breath

Pauling (1971) - used GC to detect 250 compounds in breath

Phillips (1999) - used GC/MS to detect 3000 compounds in breath

2000 - present - new advances in breath analysis each year through laser spectroscopy, mass spectrometry **and eNose analysis**

Select Publications of the Handheld eNose 2002 - 2012

University of Pennsylvania uses handheld eNose for pneumonia, CSF and sinusitis research in the ER and outpatient clinic

Lai, S. et al. Identification of upper respiratory bacterial pathogens with the electronic nose. Laryngoscope. 112(6) 975-9 (2002); Aronzon, A. et al. Differentiation between cerebrospinal fluid and serum with electronic nose. Otolaryngol Head Neck Surg 133(1) 9-16 (2005); Thaler, ER, Hanson, CW. Use of an electronic nose to diagnose bacterial sinusitis. American Journal of Rhinology, 20(2) 170-172 (2006)

Cleveland Clinic uses handheld eNose for lung cancer research

Machado, R. et al. Detection of Lung Cancer by Sensor Array Analyses of Exhaled Breath. Am. J. Respir. Crit. Care Med. 171(11) 1286-1291 (2005)

University of Amsterdam uses handheld eNose for COPD, asthma and cancer


Dragonieri, S. et al. An electronic nose in the discrimination of patients with asthma and controls. J Allergy Clin Immunol. 120(4): 856-62 (2007); Dragonieri, S. et al. An electronic nose in the discrimination of patients with non-small cell lung cancer and COPD. Lung Cancer 64(2) 166-70 (2009); Lazar, Z. et al. Electronic nose breathprints are independent of acute changes in airway caliber in asthma. Sensors 10(10) 9127-38 (2010); Fens, N. et al. Exhaled air molecular profiling in relation to inflammatory subtype and activity in COPD. Eur Respir. J. 38 1301-1309 (2011).

Research teams around the world use the handheld eNose: *Australia, New Zealand, Germany, Hungary, Italy ...*

Dent, A. et al. Electronic nose distinguishes lung cancer from healthy smoking controls. Supplement to Journal of Thoracic Oncology: (2010); Hattesoehl AD, et al. Discrimination between COPD patients with and without alpha 1-antitrypsin deficiency using an electronic nose. Respiratory Disease Diagnosis 16(8) 1258-64 (2011); Chapman, EA et al. A breath test for malignant mesothelioma using an electronic nose. Eur Respir J. December 2011. Timms, C. et al. Detection of gastro-oesophageal reflux disease (GORD) in patients with obstructive lung disease using exhaled breath profiling. J Breath Res. January 2012.

Breath Biomarkers


Disease	Compound as a disease marker	Analysis Instrument
Acute cardiac allograft rejection	Pentane	GC/FID
Myocardial infarction (MI)	Hydrocarbons	GC/FID
Asthma	Nitric Oxide	CL analyzer
COPD / ARDS	NO, CO	CL analyzer
Breast Cancer	Pentane	GC/FID
Diabetes	Acetone	GC/FID
Hemolysis	Carbon monoxide	EC CO analyzer GC/TCD
H. pylori infection	$^{13}\text{CO}_2$ or $^{14}\text{CO}_2$	Isotope Ratio MS Isotope Ratio IR
Alcoholic liver disease	Pentane	GC/FID
Liver cirrhosis	Dimethyl sulfide Volatile fatty acid	GC/FPD GC/FID
Weight Reduction	Acetone	GC/FID



VOCs from bacterial infection
acids, alcohols, aldehydes, amines, ketones,
hydrocarbons, sulfur compounds

Example:


Free radicals produce
measurable volatile
products of
oxidative stress



semivolatiles in
breath condensate



volatiles excreted
in breath

Breath Chemical Analysis



GC/MS can identify many but not all breath constituents

Breath Collection




air supply

1L sample can

Elaborate means for collection of breath constituents for analysis

Why use an “Electronic Nose”?



- Measure the entire breath profile
- Simplicity of sample measurement
- An eNose received FDA approval in 2002 for
Urinary Tract Infection (UTI)
Bacterial Vaginosis (BV)

Selected Medical Applications

▶ Bacteria Identification

- ENT bacteria, infant Otitis Media, adult Urinary Tract Infection

▶ Univ. Pennsylvania Hospital (HUP)

- Ventilator Associated Pneumonia
- Sinusitis
- Cerebrospinal fluid in the ER

▶ Cleveland Clinic Foundation (CCF)


- Lung Cancer
- ARDS, COPD, asthma
- CF, PPH

▶ Univ. of Amsterdam (AMC)

- Asthma
- Small cell cancer and COPD

Identification of ENT Bacteria



Canonical Discriminant Analysis



2 Week Prediction Success


	tests	ID
strep A+B	28/28	100%
Staph	13/13	100%
H. flu	15/15	100%
B. catarrh.	15/15	100%

Identification of ENT Bacteria



Lai SY et al, Identification of Upper Respiratory Bacterial Pathogens with the Electronic Nose. *Laryngoscope* 112, 975-9 (2002).

Identification of Eye Bacteria



Bacteria classification using Cyranose 320 electronic nose. Biomed Eng Online. 2002; 1 (1): 4
Ritaban Dutta,¹ Evor L. Hines,¹ Julian W. Gardner,¹ and Pascal Boilot¹

Breath Signature of Bacterial Infection

VOCs from bacterial metabolism

alcohols
ketones
aldehydes

organic acids
hydrocarbons
sulfides, thiols
amines



O₂

amino acids
carbohydrates

Compounds from immune system response to infection


Inflammation:

NO_x CO VOCs

Enterotoxin stimulus:


VOCs ?

Bacteria	Metabolites
<i>S. aureus</i>	2-methylbutanol 3-methylbutanol
<i>Strept. pneumoniae</i>	2-butanol lactic acid
<i>H. influenzae</i>	acetic acid indole



Ventilator Associated Pneumonia (VAP)

Cyranose used to measure exhaled breath in ICU at UPenn Hospital



Cyranose
320




Sample
taken from
exhaled
gases.


Ventilator Associated Pneumonia



William Hanson, MD
Erica Thaler, MD



eNose measurements on exhaled breath compare favorably to the combined pulmonary infection score (CPIS) used to confirm ventilator associated pneumonia in the ICU



Amer. Thoracic Society 2002



Acute Rhinosinusitis



Model	c=100, w=0.5		c=10, w=5	
	# correct	% correct	# correct	% correct
SVM	123/123	100	118/123	95.9
SVM+PCA(2)	123/123	100	113/123	91.9
SVM+PCA(3)	123/123	100	121/123	98.4

Nasal breathing cup

22 subjects

11 neg. controls

11 positives

4 months

Methods

- Nasal swabs
 - sampling of infection hotspot with calgiswab
- Nasal breathing cup

Respiratory Disease




Cleveland Clinic Foundation

Lung Cancer
Asthma
COPD
ARDS
CF

Serpil Erzurum, MD
Raed Dweik, MD
Roberto Machado, MD



Lung Cancer




Machado RF et al (CCF) Detection of Lung Cancer by Sensor Array Analyses of Exhaled Breath. *Amer. Jour. Respiratory & Critical Care Medicine* v171 1286-1291 (2005).

- Detection of lung cancer in non-smokers is feasible
- Discrimination from several disease controls: COPD, ARDS, PPH, asthma, a-1, CBE
- Goal: early detection of small tumors

DiNatale, et al (Univ Rome) Lung cancer identification by the analysis of breath by means of an array of non-selective gas sensors. *Biosensors & Bioelectronics* v18 1209-1218 (2003)

Cystic Fibrosis



- Discrimination of chronically colonized CF subjects is feasible
- In preliminary tests, 93.4% of breath samples were identified correctly

Conclusions & Hypotheses

- ▶ Many diseases produce a measurable pattern of volatile chemicals in breath, urine and blood
- ▶ **Non-invasive** breath measurement will provide rapid **diagnosis** and **treatment monitoring** capability for physicians in emergency and point-of-care applications
- ▶ **Low cost** and **low power** intelligent sensor array devices will enable home health diagnosis and monitoring capability for many individuals



For application or product information,
please contact:

Sensigent

sales@sensigent.com

626-768-2626

www.sensigent.com

